



Endo Lithotripter

Used for Purpose Of C.B.d. Stone Retrieval through the Scope.
And for Lithotripsy to Crush Stones in Biliary Duct.

Features-

Handle

Locking Feature

S.S. Sheath & Tip

S.S. Sheath Locking to Handle PFE Sheath

- Light weight with easy grip and smooth rolling of Knob.
- Unidirectional movement locking system to facilitate secure crushing / retrieval Of big C.B.D. Stone.
- Inner cutting edges to perform lithotripsy .Outer tip rounded a traumatic. Outer S.S.sheath offering no strain on Duodenoscope elevation channel.
- Secure and firm fitting of S.S. Sheath with handle
- Covers Basket and enters in C.B.Duct through Instrument channel with low Friction
Easy glide from outer s.s sheath outer S.S.sheath offering no strain on Duodenoscope elevation channel.



Technical Information-

- Working Length- 210 cms.
- Compatible with instrument Channel- 4.2 mm.
- Basket Diameter- 35 mm
- Basket Shape-Hexagonal
- Flush Port for Contrast- Flush.
- Atraumatic basket tip for ease of insertion and procedure.

Advantage-

No need to remove endoscope for lithotripsy Through the Scope-under vision C.B.D.Stone lithotripsy-Retrieval procedure done.

Availability- The Device is supplied Non- Sterile – Reusable. Flexible Sheath and Filament with Basket is available separately.

Ordering Information- - Mention product code x quantity required

Product Name	Product Code	Sheath dia.	Sheath Length cms
Shaili's Endo Lithotripter	S E L	3.7 MM	210 CMS.
Shaili Endolithotripter Flexible S.S. Sheath only	S –EL –SHEATH	3.7 MM	
Shaili Endolithotripter Basket with filament only	S- EL – BASKET		

Contact-

E-mail - drparesh@shailitraders.com www.shailitraders.com,
Telephone- +91-265 -2353153 Fax- +91-265-2350776 ,

Distributed By :

