



ISO:13485:2003
ISO:9001:2000



Foreign Body Holder Mono Filament Memory Basket

Used to retrieve foreign body from G.I.Tract - endoscopically.

Technical Information-

- Three basket sizes to facilitate Foreign Body Retrieval from G.I.Tract.
- Basket shapes - Oval, Hexagonal shaped basket.
- Retains shape for longer time.
- Handle - Light weight with Flush port for Contrast - Drug - Flush.
- Not Recommended for Lithotripsy.

Availability- Non Sterile - Reusable.



Ordering Information- Mention the product code x quantity required.

Product Name	Product Code	Basket Diameter Ø	Compatible Inst. Channel Diameter Ø	Working length cms.
Monofilament Memory Basket - Big Oval	SMMB -B	45 mm.	2.8 mm.	230 cms.
Monofilament Memory Basket - Medium Oval	SMMB -M	35 mm.	2.8 mm.	230 cms.
Monofilament Memory Basket - Hexagonal	SMMB -H	35 mm.	2.8 mm.	230 cms.

Contact-

E-mail - drparesh@shailitraders.com www.
shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :

Foreign Body Holder
Mono Filament Memory Basket

