




## Biliary Dilator

Used to dilate papilla or biliary strictures endoscopically.

### Technical Information-

- Tapering length- 3 cms.
- Tip Diameter Ø - 4fr. / 6 fr.
- SIZE- Biliary Dilator 7 fr. / 10 Fr.
- Radiopaque PTFE material .
- Availability.
- Sterile - Single use only.

**Individually packed.**



**Ordering Information-** Mention the product code x quantity required.

Product Name	Product Code	Tip Diameter Ø	Compatible with Instrument Channel mm. Diameter Ø	Working length cms.
Biliary Dilator 7 Fr.	S - BD- 7	4 Fr.	3.2 mm instrument channel	200 cms.
Biliary Dilator 10 Fr.	S - BD- 10	6 Fr.	4.2 mm instrument channel	200 cms.

### Contact-

**E-mail - [drparesh@shailtraders.com](mailto:drparesh@shailtraders.com)** www.shailtraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

**Distributed By :**

Biliary Dilator

