




# Biliary Stents

Used to drain bile from obstructed biliary duct.

## Technical Information-

- 1] Diameter  $\varnothing$  7Fr and 10 Fr. Size.
- 2] Length 5 cm. to 15 cm. as per customer requirement.
- 3] Material choice- Radiopaque PTFE AND PU.
- 4] Varients- Straight S  
Double Pig Tail DP  
Amsterdam AMST  
Standard STD
- 5] Available - Sterile, Single use.
- 6] Wireguide Compatibility - 0.038".



**Ordering Information-** Mention the product code x required Length x QTY required.

### Specifications

Material-

Product Code

### Varients-

1- Straight

2-Double Pig Tail

3-Standard

4-Amsterdam

### Biliary Stent 7 Fr.

PTFE

PU

SBS-TEF-7

SBS-PU-7

SBS-TEF-7-S

SBS-PU-7-S

SBS-TEF-7-DP

SBS-PU-7-DP

SBS-TEF-7-STD

SBS-PU-7-STD

SBS-TEF-7-AMST

SBS-PU-7-AMST

### Biliary stent 10Fr.

PTFE

PU

SBS-TEF-10

SBS-PU-10

SBS-TEF-10-S

SBS-PU-10-S

SBS-TEF-10-DP

SBS-PU-10-DP

SBS-TEF-10-STD

SBS-PU-10-STD

SBS-TEF-10-AMST

SBS-PU-10-AMST

### Contact-

**E-mail - [drparesh@shailtraders.com](mailto:drparesh@shailtraders.com)** www.  
shailtraders.com, Telephone- +91-265-2353153 Fax-+91-265-

**Distributed By :**

Biliary Stents