




Self Expandable Metallic Biliary Stent

Self Expandable Metallic Biliary Stent

Used for Stenosis or Obstructed Biliary tract caused by benign or malignant changes.

- With Introducer System.
- Radio therapy compatible.
- Stent Diameter. ϕ - 10mm.
- Stent Wire Diameter. ϕ - 0.16mm.
- Stent Length - 3 Sizes - 60mm - 80mm - 100mm.
- Uncovered stent.
- Stent introducer system Diameter. ϕ - 8 fr.
- Compatible with 0.038 inch OD ϕ wire guide.
- Compatible with instrument channel size - 3.2mm instrument channel Diameter. ϕ .
- Availability - Individually packed. Sterile - For single use only.



Ordering Information- Mention the product code x quantity required

Product Name	Product Code	Stent Length mm	Stent Diameter ϕ mm
Self Expandable Metallic Biliary Stent	S - MBS - 6	60mm.	10mm.
Self Expandable Metallic Biliary Stent	S - MBS - 8	80mm.	10mm.
Self Expandable Metallic Biliary Stent	S - MBS - 10	100mm.	10mm.

Contact-

E-mail - drparesh@shailitraders.com www.shailitraders.com, Telephone- +91-265 -2353153 Fax- +91-265-

Distributed By :

