




Biopsy Forceps

Used Endoscopically to obtain Mucosal Tissue Samples from Gastrointestinal Tract for Microscopic - Histopathological Examination.

Technical Information-

- Type - With Spike / Without Spike.
- Availability- Sterile, For Single use only.
- Compatible with 2.8 mm. instrument channel.
- Length- 160 cms , 210 cms.
- Individually packed.



Ordering Information- Mention the product code x quantity required.

| Product Name | Product Code | Type | Compatible with Instrument channel mm. Diameter Ø | Working length cms. |
|----------------|--------------|----------------|---------------------------------------------------|---------------------|
| Biopsy Forceps | S- BF- WS-16 | With Spike | 2.8 mm. | 160 cms. |
| Biopsy Forceps | S- BF- WS-21 | With Spike | 2.8 mm. | 210 cms. |
| Biopsy Forceps | S- BF- -16 | With Out Spike | 2.8 mm. | 160 cms. |
| Biopsy Forceps | S- BF- -21 | With Out Spike | 2.8 mm. | 210 cms. |

Contact-

E-mail - drparesh@shailitraders.com www.shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :



Biopsy Forceps