




Cleaning Brush

Used to clean endoscope channel .

Technical Information-

- Type - Both end bristles / One end bristles.
- Shaft - Plastic coated wire / S.S. Coil.
- Semi - Soft and channel friendly bristles for non invasive cleaning of endoscope channel.
- Tip - Atraumatic tip.
- Availability- Non - Sterile.
- Single use / Reusable.
- Compatible with 2.8 mm. instrument channel.
- Length- 200 cms.
- Individually packed.



Ordering Information- Mention the product code x quantity required.

Product Name	Product Code	Type	Compatible with Instrument channel mm Diameter Ø	Working length cms.
Cleaning Brush	S- CB-200	Both end bristles.	2.8 mm.	200 cms.
Cleaning Brush	S- CB-S-200	One end bristles	2.8 mm.	200 cms.

Contact-

E-mail - drparesh@shailitraders.com

www.shailitraders.com, Telephone- +91-265 -2353153

Fax-+91-265-2350776 ,

Distributed By :

Cleaning Brush

