



# ERCP Catheter/Canula

**Used endoscopically to opacify the biliary duct with radiopaque contrast for fluoroscopic view.**

## Technical Information-

Type - Tapered tip - Ball tip

Catheter tip Diameter  $\emptyset$  - 5 Fr.

Recommended wire guide Diameter  $\emptyset$  0.035 inch Diameter  $\emptyset$  Flexible.

Availability - Sterile For single use only.

Length- 200 cms.

**Individually packed. .**



**Ordering Information-** Mention the product code x quantity required.

Product Name	Product Code	Type	Compatible with Instrument channel mm Diameter $\emptyset$	Working length cms.
ERCP CATHETER	S ERC-T	TAPERED TIP	3.2 mm.	200 cms.
ERCP CATHETER	S ERC-B	BALL TIP	3.2 mm.	200 cms.

## Contact-

**E-mail - [drparesh@shailitraders.com](mailto:drparesh@shailitraders.com)** www.shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

**Distributed By :**

Ercep Catheter

