




ISO:13485:2003  
ISO:9001:2000



## Esophageal Dilators-

Used to Dilate Esophageal Strictures over previously placed wire guide.

- Non Sterile and Reusable.
- Available Sizes - 5mm to 16mm Outer Dia.  $\phi$



- Soft, Flexible and Radiopaque Tip & Radiopaque Junction.
- Gradual Tapering.
- Low friction lumen surface.
- Compatible with 0.038 inch OD.  $\phi$  wire guide.



### Ordering Instructions -

**Dilator size** - Diameter - 5 mm. To 16 mm. Length Two sizes available -100 cms , 70 cms.

**Availability** - 1] Set of 6 Pcs. Even Dia.  $\phi$  Size . / Odd Dia.  $\phi$  Size 2] Set of 12 Pcs. 3] Individual Pcs. Of desired Dia.  $\phi$   
Set includes 1x Shaili Flexible spring tip Guide Wire 200 cms long. Carrying case is available with order of set only.

### Esophageal Dilators-

Product Code	Outer Diameter $\phi$	Length	Product Code	Outer Diameter $\phi$	Length
SED-500	5 mm.	100 cms.	SED-570	5 mm.	70 cms.
SED-600	6 mm.	100 cms.	SED-670	6 mm.	70 cms.
SED-700	7 mm.	100 cms.	SED- 770	7 mm.	70 cms.
SED-800	8 mm.	100 cms.	SED-870	8 mm.	70 cms.
SED-900	9 mm.	100 cms.	SED-970	9 mm.	70 cms.
SED-1000	10 mm.	100 cms.	SED-1070	10 mm.	70 cms.
SED-1100	11 mm.	100 cms.	SED-1170	11 mm.	70 cms.
SED-1200	12 mm.	100 cms.	SED-1270	12 mm.	70 cms.
SED-1300	13 mm.	100 cms.	SED-1370	13 mm.	70 cms.
SED-1400	14 mm.	100 cms.	SED-1470	14 mm.	70 cms.
SED-1500	15mm.	100 cms.	SED-1570	15mm.	70 cms.
SED-1600	16 mm.	100 cms.	SED-1670	16 mm.	70 cms.

**Ordering Information-** Mention the product code x quantity required.

**Distributed By :**

**Contact-**

**E-mail - [drparesh@shailitraders.com](mailto:drparesh@shailitraders.com) www.**

shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Esophageal Dilators

