



## Esophageal Self Expandable Metallic Stent

Used to treat Tracheoesophageal Fistulas and expansion treatment of the Esophageal Stenosis Caused by Esophageal Carcinoma, Chemical Injuries, Cardiac Carcinoma, Anastomatic stenosis. And Other Traumas, which are surgical contraindicated.

### Technical Information

- Radio therapy compatible.
- Esophageal Stents with Introducer system.
- Stent Diameter.  $\phi$  - 18mm.
- Stent Wire Diameter.  $\phi$  - 20mm.
- Stent Length - 4 Sizes - 50mm - 100mm - 120mm - 150mm covered stent.
- Compatible with 0.038 inch OD  $\phi$  wire guide.
- Availability - Individually packed. Sterile - For single use only.



**Ordering Information-** Mention the product code x quantity required.

### Product Name

Esophageal Self Expandable Metallic Stent Covered  
Esophageal Self Expandable Metallic Stent Covered  
Esophageal Self Expandable Metallic Stent Covered  
Esophageal Self Expandable Metallic Stent Covered

### Product Code

S - ES - 5  
S - ES - 10  
S - ES - 12  
S - ES - 15

### Stent Length mm

50mm  
100mm  
120mm  
150mm

### Stent Diameter mm

18mm.  
18mm.  
18mm.  
18mm.

### Contact-

E-mail - [drparesh@shailitraders.com](mailto:drparesh@shailitraders.com) www.  
shailitraders.com, Telephone- +91-265-2353153 Fax- +91-265-

### Distributed By :



Esophageal Self Expandable Metallic Stent

