



## Flexible Spring Tip Metalic Guide Wire

**Used to introduce Esophageal Dilators Endoscopically for Dilatation.**

- Specially designed Flexible atraumatic tip Diameter  $\emptyset$  - 2.0 mm.
- Spring base and metallic wire shaft are interlocked.
- Joint of Spring and wire shaft is tapered to stop guide wire stacking in the dilator.
- Smooth surface finish offers easy passage through dilator lumen during dilatation procedure.

- **Availability- Non Sterile Reusable.**
- **Individually packed in plastic protector sheath.**



**Ordering Information-** Mention the product code x quantity required.

Product Name	Product Code	Tip Diameter $\emptyset$	Compatible with Esophageal dilator inner lumen Diameter $\emptyset$	Working length cms.
Flexible Spring Tip Guide Wire	SMGW	2.0 mm.	1.0 mm.	200 cms.

**Contact-**

**E-mail - [drparesh@shailitraders.com](mailto:drparesh@shailitraders.com)** www.  
shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

**Distributed By :**

Flexible Spring Tip Marked  
Metalic Guide Wire

