



ISO:13485:2003
ISO:9001:2000



Foreign Body Holding Forcep 2 Prong, 3 Prong And Loop Basket

Used to retrieve Foreign body from G.I. Tract.

Technical Information-

Compatible with 2.8 mm. instrument channel. Working length 170 cms.

Supplied Non Sterile and are Reusable.

- Individually packed in plastic protector sheath . Excellent and secure grip to retrieve Foreign Body from G.I. Tract.
- 2 prong and 3 prong and Loop Basket .

* Available for 1.8 mm instrument channel diameter against confirm order.



Ordering Information- Mention the product code x quantity required.

Product Name	Product Code	Working Lengthcms	Compatible with Instrument Channel mm. Diameter Ø	Jaws
F.B.Holder 2 prong	SFBH- 223	170 cms.	2.8 mm. Dia.	2 prong
F.B.Holder 3 prong	SFBH- 323	170 cms.	2.8 mm. Dia.	3 prong
F.B.Holder 2 prong	SFBH- 218	170 cms.	1.8 mm. Dia.	2 prong
F.B.Holder 3 prong	SFBH- 318	170 cms.	1.8 mm. Dia.	3 prong
F.B.Holder Loop Basket	SFBH- LB	170 cms..	2.8 mm. Dia.	Loop Basket


Shaili Range of Foreign Body Holders -

2 prong, 3 Prong , Loop Basket , Retrieval Basket , Snare Shark Tooth Forcep, Alligator Forcep
Monofilament Memory Basket.

Contact-

E-mail - drparesh@shailitraders.com www.shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :



Foreign Body Holding Forcep

