



ISO:13485:2003  
ISO:9001:2000



Foreign Body Holding Forceps

## Foreign Body Holding Forceps

Shark Tooth Forcep.  
Alligator Forcep.

### Features-

- Used to retrieve foreign body from G.I. Tract.
- Can be used in ERCP procedures as well.
- Individually packed in plastic protector sheath .
- Supplied Non Sterile and it is Reusable.
- Specially produced outer sheath for easy elevation through elevator.
- Firm grip & easy movement.



**Ordering Information-** Mention the product code x quantity required.

No.	Product Code	Name	Outer Diameter Ø	Inst.Channel Diameter Ø	Jaws Specifications Working Length cms.
1	SSTF-230	SHARK TOOTH FORCEPS	2.5 mm.	2.8 mm.	Rat Tooth Jaws 230 cm.
2	SATF-170	ALLIGATOR FORCEPS	2.5 mm.	2.8 mm.	.Aligator Jaws 170 cm.

**Shaili Range of Foreign Body Holder :** 2 prong , 3 Prong , Retrieval Basket , Snare, Shark Tooth Forcep  
Alligator Forcep, Monofilament Memory Basket, Loop Basket.

### Contact-

**E-mail - [drparesh@shailitraders.com](mailto:drparesh@shailitraders.com)** www.  
shailitraders.com, Telephone- +91-265 -2353153 Fax- +91-265-

**Distributed By :**

