



Manual Lithotripter

Used For C.B.D. Stone Retrieval Lithotripsy.

Technical Information-


Handle- Light weight with unidirectional movement.

Sheath tip- Inner cutting edges to perform lithotripsy. Outer tip rounded atraumatic.

Secure and firm fitting of S.S. Sheath with handle to facilitate retrieval - lithotripsy procedure.

Availability- Non Sterile.

Reusable.



Ordering Information- Mention the product code x quantity required.

Product Name	Product Code	Sheath diameter Ø	Sheath Length cms.
Shaili's Manual Lithotripter	SML	4.7 mm.	95 cms.

Contact-

E-mail - drparesh@shailitraders.com www.shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :

Manual Lithotripter

