





Needle Knife

Used for Precut Papillotomy.

Technical Information-

- Tapered tip.
- Working length 200 cms.
- Cathetar diameter \emptyset - 5 fr.
- Individually packed.
- Supplied Sterile and it is for single use.
- Active Cord available separately.



Ordering Information- Mention the product code x quantity required.

No.	Product code	Length cms	Needle length	Outer Dia. \emptyset	Cutting wire
1	SNK-100	200 cms.	5 mm.	1.70 mm.	S.S.Mono Filament.

Contact-

E-mail - drparesh@shailtraders.com www.
shailtraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :

Needle Knife

