



Neg Wire Guide

**Used to gain access through difficult biliary and pancreatic strictures .
To facilitate exchange of compatible wire guided accessories to the
desired biliary or pancreatic ducts.**

Technical Information-


Size - Three OD. Ø Sizes 0.032 inch , 0.035 inch. , 0.038 inch.

Availability - Individually packed. Sterile For single use only.

Compatible with 2.8 mm instrument channel.

Length - 400 cms.

Tip - Flexible tip .



Ordering Information- Mention the product code x quantity required.

Product Name	Product Code	Compatible with Instrument channel mm. Diameter Ø	Working length cms.
NEG WIRE GUIDE - 0.032 inch.	S -NEG - 032	2.8 mm.	400 cms.
NEG WIRE GUIDE - 0.035 inch.	S -NEG - 035	2.8 mm.	400 cms.
NEG WIRE GUIDE - 0.038 inch.	S -NEG - 038	2.8 mm.	400 cms.

Contact-

E-mail - drparesh@shailtraders.com www.
shailtraders.com, Telephone- +91-265 -2353153 Fax- +91-265-

Distributed By :

Neg Wire Guide

