





Polypectomy Snare

Used For Endoscopic Excision Of Polyps Using Monopolar Electro-surgical Unit And Retrieval Of Foreign Bodies From G.I. Tract.

Technical Information-

- Availability - Single use only.
- With flush port to clean wire.
- Outer sheath - PTFE.
- Shape - Hexagonal / Oval.

Individually packed .



Ordering Information- Mention the product code x quantity required.

Product Name	Product Code	Loop Shape	Compatible With Instrument Channel mm. Diameter Ø	Working length cms.
Polypectomy snare	SPSO	OVAL	2.8 mm.	220 cms.
Polypectomy Snare	SPSH	Hexagonal	2.8 mm.	220 cms.

Contact-

E-mail - drparesh@shailtraders.com www.shailtraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :



Polypectomy Snare