



ISO:13485:2003
ISO:9001:2000





Foreign Body Holding Forcep Retrieval Basket

Foreign Body Holding Forcep Retrieval Basket

Used to retrieve Foreign Bodies and C.B.D. stones.

Technical Information-

- Working length 230 cms.
- Compatible with instrument channel 2.8 mm. / 3.2 mm.
- Basket Diameter- Available in 3 sizes - 30 mm. - 35 mm. - 40 mm.
- Multifilament, Hexagonal shaped.
- Flush port for contrast-flush.
- Atraumatic Basket Tip for ease of insertion and procedure.
- User friendly handle for Lithotripsy.
- Compatible with Lithotripter.
- Reusable, Non Sterile.
- Individually packed.



Ordering Information- For 3.2mm. or above channel size, Strand diameter \emptyset is 55mm.
For 2.8mm. Channel Size, strand diameter \emptyset is .41mm.
Working length 2300mm.

Product Name	Product Code	Basket Diameter \emptyset	Instrument Channel Size \emptyset mm.	Sheath Diameter \emptyset
Shaili Retrieval Basket	S R B - 4302	30 mm.	2.8 mm.	2.50mm.
Shaili Retrieval Basket	S R B - 4352	35 mm.	2.8 mm.	2.50mm.
Shaili Retrieval Basket	S R B - 4402	40 mm.	2.8 mm.	2.50mm.
Shaili Retrieval Basket	S R B - 4303	30 mm.	3.2 mm.	2.50mm.
Shaili Retrieval Basket	S R B - 4353	35 mm.	3.2 mm.	3.00mm.
Shaili Retrieval Basket	S R B - 4403	40 mm.	3.2 mm.	3.00mm.

Shaili Range of Foreign Body Holders- 2 prong , 3 Prong , Retrieval Basket , Snare , Shark Tooth Forcep Alligator Forcep , Loop Basket, Monofilament Memory Basket.

Contact-

E-mail - drparesh@shailitraders.com www.
shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :