



Special Purpose Injector Needle

Injector Needles- Sclero Therapy Needle

Used for Endoscopic Injection in to Esophageal Varices above O.G.Junction.

Technical specifications -

Compatible with Instrument channel Needle	2.8 mm. instrument channel.
Working length cms.	Short and Sharp Bevel, Available in 21 swg. & 24 swg. 170 cms.
Availability	Sterile- Individually packed .Single use only.

- Specially developed needle secure firm grip with inner tube & possible leakage from the joint of inner tube & needle.

- Specially developed needle give "JET" pressure to inject the drug.



Ordering Information- Mention the product code x quantity required.

No.	Product code	Product name	Needle gauge	LengthCms	Instrument Chanel Diameter Ø
01	SGN - 2321	Sclerotherapy Needle	21 G	170 cms.	2.8 mm.
02	SGN- 2321- C	Sclerotherapy Needle	21 G	200 cms.	2.8 mm.
03	SSN - 2324	Sclerotherapy Needle	24 G	170 cms.	2.8 mm.
04	SSN- 2324-C	Sclero therapy Needle	24 G	200 cms.	2.8 mm.

Contact-

E-mail - drparesh@shailitraders.com www.shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :

Special Purpose Injector Needle

