




# Sphincterotome

Used for Endoscopic Cannulation of Biliary Ductal System and Sphincterotomy.

## Technical Information-

- Triple Lumen with Precurved tapered tip.
- Working length 200 cms.
- Cutting wire with **protective insulation sheath**, unique protective feature for endoscope.
- **Unique guide wire port to pass the wire straight without any curve & less friction on guide wire coating.**
- Supplied Sterile and it is for single use.
- Active Cord available separately.



**Ordering Information-** Mention the product code x quantity required.

	Product code	Tip length	Working Length	Insulated Wire	Cutting Wire	Wire loop	Wire guide compatibility
1	SSPH-SN	5mm.	200 cms.	8mm.	S. S.Monofilament	25mm.	0.038"
2	SSPH -LN	10mm.	200 cms.	8mm.	S.S.Monofilament	25mm.	0.038"
3	SSPH -NN	0mm.	200 cms.	8mm.	S.S.Monofilament	25mm.	0.038"

## Contact-

E-mail - [drparesh@shailtraders.com](mailto:drparesh@shailtraders.com) www.  
shailtraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :

Sphincterotome

