



ISO:13485:2003
ISO:9001:2000



- **Band Ligator Set.**
- **Barrel With 6 Loaded Bands.**

Used for Endoscopic Ligation of Esophageal Varices at or above Gastro Esophageal Junction.


This device is supplied Non Sterile, Individually packed - Intended for Single use only.

Technical Information-

Shaili Band Ligator requires 9.5 to 11.2 .mm endoscope diameter Ø.

Band Ligator Set Contains-

Band Ligator x1 pc., Barrel loaded with 6 x Bands, Atraumatic plastic hook with cannula. X1 pc., Flushing Needle blunt tip x1 pc.



Ordering Information- Mention the product code x quantity required.

Product

Band Ligator Set
Barrel Only with 6 loaded bands

Product Code

SBLS-6
SBL-6

Contact-

E-mail - drparesh@shailitraders.com www.shailitraders.com, Telephone- +91-265 -2353153 Fax-+91-265-

Distributed By :

Band Ligator Set
Barrel With 6 Loaded Bands

