



# Wire Guided Stone Extraction Balloon Catheter

**Used For Endoscopic Extraction Of Biliary Stones .**

## Technical Information-

**Double lumen feature** provides wire guide access to facilitate stone extraction from desired location during procedure.

Catheter Diameter  $\emptyset$  - 7 Fr.


Compatible with 0.038 inch OD wire guide.

Availability - Individually packed. Sterile For single use only.

Working Length - 200 cms.

Balloon Diameter  $\emptyset$  - 14 mm.

Balloon Length - 16 mm.



**Ordering Information-** Mention the product code x quantity required.

Product Name	Product Code	Baloon length mm.	Baloon Diameter $\emptyset$	Max Air $\text{Cm}^2$	Max. Liquidml	Catheter Diameter $\emptyset$
Double Lumen Stone Extraction Baloon Catheter	S- SEBC - 7	16 mm	14 mm.	4.5 $\text{Cm}^2$	2.5 ml.	7 Fr.

## Contact-

**E-mail - [drparesh@shailtraders.com](mailto:drparesh@shailtraders.com)** www.shailtraders.com, Telephone- +91-265-2353153 Fax-+91-265-

**Distributed By :**

Wire Guided  
Stone Extraction Balloon Catheter

